

For Research Use Only

CNR1 Recombinant antibody, PBS Only (Detector)

Catalog Number: 82787-3-PBS



Basic Information

Catalog Number:

82787-3-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG12625

GenBank Accession Number:

BC074812

GeneID (NCBI):

1268

UNIPROT ID:

P21554

Full Name:

cannabinoid receptor 1 (brain)

Calculated MW:

472 aa, 53 kDa

Purification Method:

Protein A purification

CloneNo.:

250285B1

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP02192-2, Cannabinoid receptor 1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82787-5-PBS. Detection antibody: 82787-3-PBS. Standard: Ag12625. Range: 1.25-20 ng/mL