

For Research Use Only

HNRNPC Recombinant antibody, PBS Only (Capture)

Catalog Number: 82798-5-PBS



Basic Information

Catalog Number:

82798-5-PBS

Size:

1mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG2356

GenBank Accession Number:

BC003394

GeneID (NCBI):

3183

UNIPROT ID:

P07910

Full Name:

heterogeneous nuclear
ribonucleoprotein C (C1/C2)

Calculated MW:

32 kDa

Purification Method:

Protein A purification

CloneNo.:

230373F6

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00088-1, HNRNPC Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82798-5-PBS. Detection antibody: 82798-4-PBS. Standard: Ag2356. Range: 1.25-20 ng/mL.



Biolayer interferometry (BLI) kinetic assay of 82798-5-PBS against Human HNRNPC was performed. The affinity constant is 0.77 nM.