

For Research Use Only

OLIG2 Recombinant antibody, PBS Only (Capture)

Catalog Number: 82806-5-PBS



Basic Information

Catalog Number:

82806-5-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG5089

GenBank Accession Number:

BC047511

GeneID (NCBI):

10215

UNIPROT ID:

Q13516

Full Name:

oligodendrocyte lineage transcription

factor 2

Calculated MW:

32 kDa

Purification Method:

Protein A purification

CloneNo.:

240839D6

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00809-1, OLIG2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82806-5-PBS. Detection antibody: 82806-6-PBS. Standard: Ag5089. Range: 0.625-80 ng/mL