

For Research Use Only

Haptoglobin Recombinant monoclonal antibody, PBS Only

Catalog Number: 82812-11-PBS



Basic Information

Catalog Number: 82812-11-PBS	GenBank Accession Number: BC121125	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 3240	CloneNo.: 25217683
Isotype: IgG	UNIPROT ID: P00738	
Immunogen Catalog Number: EG6628	Full Name: haptoglobin	
	Observed MW: 42 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
human

Background Information

HP(Haptoglobin) is also named as zonulin and belongs to the peptidase S1 family. HP, a plasma glycoprotein that binds free hemoglobin, has a tetrameric structure of 2 alpha(16 kDa and 9 kDa) and 2 beta(40 kDa) polypeptides that are covalently associated by disulfide bonds. In most species, apart from ruminants, Hp has a molecular mass of 100 kDa, consisting of two subunits of 40 kDa and two subunits of 9 kDa, although in a few species, such as man, genetic variant of Hp forms polymers of higher mass(PMID:2361363). Recent studies of haptoglobin show that certain oligosaccharide structures predominate in different diseases. For example, a highly-fucosylated structure is found in breast cancer and ovarian cancer, highly-sialylated structures in Crohn's disease and highly branched structures in alcoholic liver disease and fucosylated haptoglobin is a good serum marker for pancreatic cancer.(PMID:16385567).

Storage

Storage:
Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:

T: 4006900926


E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data

Various lysates were subjected to SDS PAGE followed by western blot with 82812-11-RR (HP antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 82812-11-PBS in a different storage buffer formulation.



Biolayer interferometry (BLI) kinetic assays of 82812-11-RR against Human Haptoglobin were performed. The affinity constant is below 1 pM.