

For Research Use Only

EDNRB Recombinant antibody, PBS Only (Detector)



Catalog Number: 82818-4-PBS

Basic Information

Catalog Number:

82818-4-PBS

Size:

1mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_000115

GeneID (NCBI):

1910

UNIPROT ID:

P24530

Full Name:

endothelin receptor type B

Calculated MW:

50 kDa

Purification Method:

Protein A purification

CloneNo.:

230334A4

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve ofMP00078-1, EDNRB Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82818-6-PBS. Detection antibody: 82818-4-PBS. Standard: SY00958. Range: 0.156-10 ng/mL



Cytometric bead array standard curve ofMP00078-2, EDNRB Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82818-5-PBS. Detection antibody: 82818-4-PBS. Standard: SY00958. Range: 0.156-10 ng/mL