

For Research Use Only

# EDNRB Recombinant antibody, PBS Only (Detector)



Catalog Number: 82818-7-PBS

## Basic Information

Catalog Number:

82818-7-PBS

Size:

1mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM\_000115

GeneID (NCBI):

1910

UNIPROT ID:

P24530

Full Name:

endothelin receptor type B

Calculated MW:

50 kDa

Purification Method:

Protein A purification

CloneNo.:

230334D1

## Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve ofMP00078-3, EDNRB Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82818-3-PBS. Detection antibody: 82818-7-PBS. Standard: SY00958. Range: 0.156-10 ng/mL