

For Research Use Only

BACH1 Recombinant antibody, PBS Only (Detector)



Catalog Number: 82842-1-PBS

Basic Information

Catalog Number: 82842-1-PBS	GenBank Accession Number: BC063307	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 571	CloneNo.: 4P21
Source: Rabbit	ENSEMBL Gene ID: ENSG00000156273	
Isotype: IgG	UNIPROT ID: O14867	
Immunogen Catalog Number: AG5128	Full Name: BTB and CNC homology 1, basic leucine zipper transcription factor 1	
	Calculated MW: 82 kDa	
	Observed MW: 100-110 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP80007-2, BACH1 Hybrid Matched Antibody Pair, PBS Only. Capture antibody: 66762-2-PBS. Detection antibody: 82842-1-PBS. Standard: Ag6011. Range: 0.313-40 ng/mL