

For Research Use Only

# Mouse Snf1lk Recombinant antibody, PBS Only

Catalog Number: 82851-3-PBS



## Basic Information

Catalog Number:

82851-3-PBS

Size:

1mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG33435

GenBank Accession Number:

NM\_010831

GeneID (NCBI):

17691

UNIPROT ID:

Q60670

Full Name:

salt inducible kinase 1

Calculated MW:

85 kDa

Observed MW:

85 kDa

Purification Method:

Protein A purification

CloneNo.:

3D2

## Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, mouse

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 82851-3-RR (Snf1lk antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 82851-3-PBS in a different storage buffer formulation.



Biolayer interferometry (BLI) kinetic assays of 82851-3-RR against Mouse Snf1lk were performed. The affinity constant is 0.686 nM.