

For Research Use Only

Snf1lk Recombinant antibody

Catalog Number: 82851-4-RR



Basic Information

Catalog Number:

82851-4-RR

Size:

1000 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG33435

GenBank Accession Number:

NM_010831

GeneID (NCBI):

17691

UNIPROT ID:

Q60670

Full Name:

salt inducible kinase 1

Calculated MW:

85 kDa

Observed MW:

85 kDa

Purification Method:

Protein A purification

CloneNo.:

3A6

Recommended Dilutions:

WB 1:1000-1:6000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

mouse

Positive Controls:

WB : mouse liver tissue, mouse kidney tissue

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 82851-4-RR (Snf1lk antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours.