

For Research Use Only

BCMA/TNFRSF17 Recombinant antibody

Catalog Number: 82855-1-RR



Basic Information

Catalog Number:

82855-1-RR

Size:

300 ug/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC058291

GeneID (NCBI):

608

UNIPROT ID:

Q02223

Full Name:

tumor necrosis factor receptor
superfamily, member 17

Calculated MW:

184 aa, 20 kDa

Observed MW:

20 kDa

Purification Method:

Protein A purification

CloneNo.:

5H7

Recommended Dilutions:

WB 1:1000-1:4000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : Ramos cells, IM-9 cells

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 82855-1-RR (TNFRSF17 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.