

For Research Use Only

FBXO44 Recombinant antibody, PBS Only (Capture)



Catalog Number: 82857-8-PBS

Basic Information

Catalog Number:	82857-8-PBS	GenBank Accession Number:	BC007832	Purification Method:	Protein A purification
Size:	1mg/ml	GeneID (NCBI):	93611	CloneNo.:	230229F7
Source:	Rabbit	UNIPROT ID:	Q9H4M3		
Isotype:	IgG	Full Name:	F-box protein 44		
Immunogen Catalog Number:	AG0977	Calculated MW:	26 kDa		

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve ofMP00020-1, FBXO44 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82857-8-PBS. Detection antibody: 82857-6-PBS. Standard: Ag0977. Range: 1.25-20 ng/mL.