

For Research Use Only

Mouse GM-CSF Recombinant antibody, PBS Only (Capture)

Catalog Number: 82861-9-PBS



Basic Information

Catalog Number:

82861-9-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM-009969

GeneID (NCBI):

12981

UNIPROT ID:

P01587

Full Name:

colony stimulating factor 2
(granulocyte-macrophage)

Purification Method:

Protein A purification

CloneNo.:

230286E11

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

mouse

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00306-1, Mouse GM-CSF Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82861-9-PBS. Detection antibody: 82861-7-PBS. Standard: Eg0512. Range: 0.625-10 ng/mL



Cytometric bead array standard curve of MP00306-2, Mouse GM-CSF Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82861-9-PBS. Detection antibody: 82861-8-PBS. Standard: Eg0512. Range: 0.625-10 ng/mL