

For Research Use Only

Mouse VEGFA Recombinant antibody, PBS Only

Catalog Number: 82864-7-PBS



Basic Information

Catalog Number:

82864-7-PBS

Size:

1mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_001287056

GeneID (NCBI):

22339

Full Name:

vascular endothelial growth factor A

Purification Method:

Protein A purification

CloneNo.:

2I18

Applications

Tested Applications:

IF/ICC, Indirect ELISA

Species Specificity:

mouse

Background Information

Vegfa is a member of the PDGF/VEGF growth factor family. It is a mitogen that specifically acts on endothelial cells and has various effects, including mediating increased vascular permeability, inducing angiogenesis, vasculogenesis, endothelial cell growth, promoting cell migration, and inhibiting apoptosis. It binds to the FLT1/VEGFR1 and KDR/VEGFR2 receptors, heparan sulfate and heparin. It may play a role in increasing vascular permeability during lactation, when increased transport of molecules from the blood is required for efficient milk protein synthesis.

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunofluorescent analysis of (4% PFA) fixed NIH/3T3 cells using Vegfa antibody (82864-7-RR, Clone: 2I18) at dilution of 1:250 and CoraLite®488-Conjugated AffiniPure Goat Anti-Rabbit IgG(H+L) (SA00013-2). This data was developed using the same antibody clone with 82864-7-PBS in a different storage buffer formulation.



Biolayer interferometry (BLI) kinetic assays of 82864-7-RR against Mouse Vegfa were performed. The affinity constant is below 1 pM.