

For Research Use Only

Mouse Il12b Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number:82865-10-PBS



Basic Information

Catalog Number: 82865-10-PBS	GenBank Accession Number: NM_001303244	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 16160	CloneNo.: 252385H4
Isotype: IgG	UNIPROT ID: P43432	
Immunogen Catalog Number: EG0299	Full Name: interleukin 12b	
	Calculated MW: 38kd	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
mouse

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean IL-12/IL-23 p40 concentration was determined to be 1.3 ng/mL in Mouse serum.



Cytometric bead array standard curve of MP00413-4, MOUSE IL-12/IL-23 p40 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82865-11-PBS. Detection antibody: 82865-10-PBS. Standard: Eg0299. Range: 0.156-20 ng/mL.