

For Research Use Only

Mouse IL-12/IL-23 p40 Recombinant antibody, PBS Only (Capture)

Catalog Number: 82865-8-PBS



Basic Information

Catalog Number:

82865-8-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_001303244

GeneID (NCBI):

16160

UNIPROT ID:

P43432

Full Name:

interleukin 12b

Calculated MW:

38kd

Purification Method:

Protein A purification

CloneNo.:

240210F8

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

mouse

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00413-1, MOUSE IL-12/IL-23 p40 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82865-8-PBS. Detection antibody: 82865-6-PBS. Standard: Eg0299. Range: 0.625-80 ng/mL