

For Research Use Only

SFTPB Recombinant antibody, PBS Only



Catalog Number: 82866-16-PBS

Basic Information

Catalog Number: 82866-16-PBS	GenBank Accession Number: BC032785	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 6439	CloneNo.: 230240D4
Source: Rabbit	UNIPROT ID: P07988	
Isotype: IgG	Full Name: surfactant protein B	
Immunogen Catalog Number: AG23605	Calculated MW: 381 aa, 42 kDa	

Applications

Tested Applications:
WB, ELISA

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse lung tissue were subjected to SDS PAGE followed by western blot with 82866-16-RR (SFTPB antibody) at dilution of 1:2500 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 82866-16-PBS in a different storage buffer formulation.