

For Research Use Only

SFTPB Recombinant antibody

Catalog Number: 82866-16-RR



Basic Information

Catalog Number:

82866-16-RR

Size:

250 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG23605

GenBank Accession Number:

BC032785

GeneID (NCBI):

6439

UNIPROT ID:

P07988

Full Name:

surfactant protein B

Calculated MW:

381 aa, 42 kDa

Observed MW:

43 kDa

Purification Method:

Protein A purification

CloneNo.:

230240D4

Recommended Dilutions:

WB 1:1000-1:5000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

Human, mouse

Positive Controls:

WB : mouse lung tissue,

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse lung tissue were subjected to SDS PAGE followed by western blot with 82866-16-RR (SFTPb antibody) at dilution of 1:2500 incubated at room temperature for 1.5 hours.