

For Research Use Only

SFTPB Recombinant antibody

Catalog Number: 82866-2-RR



Basic Information

Catalog Number: 82866-2-RR	GenBank Accession Number: BC032785	Purification Method: Protein A purification
Size: 386.8 µg/ml	GeneID (NCBI): 6439	CloneNo.: 230240F4
Source: Rabbit	UNIPROT ID: P07988	
Isotype: IgG	Full Name: surfactant protein B	
Immunogen Catalog Number: AG23605	Calculated MW: 381 aa, 42 kDa	

Applications

Tested Applications:
ELISA

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data