

For Research Use Only

# SFTPB Recombinant antibody

Catalog Number: 82866-4-RR



## Basic Information

Catalog Number:

82866-4-RR

Size:

953.7 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG23605

GenBank Accession Number:

BC032785

GeneID (NCBI):

6439

UNIPROT ID:

P07988

Full Name:

surfactant protein B

Calculated MW:

381 aa, 42 kDa

Purification Method:

Protein A purification

CloneNo.:

230240E6

## Applications

Tested Applications:

ELISA

Species Specificity:

Human

## Background Information

### Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

