

For Research Use Only

NOX2 Recombinant antibody, PBS Only

Catalog Number: 82872-1-PBS



Basic Information

Catalog Number:	GenBank Accession Number:	Purification Method:
82872-1-PBS	BC032720	Protein A purification
Concentration:	GeneID (NCBI):	CloneNo.:
1mg/ml	1536	1B17
Source:	UNIPROT ID:	
Rabbit	P04839	
Isotype:	Full Name:	
IgG	cytochrome b-245, beta polypeptide	
Immunogen Catalog Number:	Calculated MW:	
AG4402	577 aa, 65 kDa	
	Observed MW:	
	55 kDa	

Applications

Tested Applications:
IF-P, FC (Intra), Indirect ELISA

Species Specificity:
human, mouse, rat

Background Information

NOX2, also named as CYBB, CGD, 91-phox, gp91-1, gp91-phox, p22 phagocyte B-cytochrome, cytochrome b-245 and beta polypeptide, is a critical component of the membrane-bound oxidase of phagocytes that generates superoxide. It is the terminal component of a respiratory chain that transfers single electrons from cytoplasmic NADPH across the plasma membrane to molecular oxygen on the exterior. This full length protein has three glycosylation sites. CYBB is found in human cardiomyocytes as multiple bands: the signal between 55 and 65 kDa is probably the unglycosylated protein, because the predicted molecular weight of unglycosylated CYBB protein is approximately 58 kDa (PMID: 17587483) to 65 kDa in phagocytes and the bands around 80 kDa probably represent glycosylated CYBB, as has also been shown in human umbilical vein endothelial cells (PMID: 12610097). In other reports, it also can be detected a band of 91 kDa (PMID: 19965781).

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:




T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunofluorescent analysis of (4% PFA) fixed paraffin-embedded mouse spleen tissue using NOX2 antibody (82872-1-RR, Clone: 1B17) at dilution of 1:2000 and CoraLite®488-Conjugated AffiniPure Goat Anti-Rabbit IgG(H+L) (SA00013-2). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 82872-1-PBS in a different storage buffer formulation.

1x10⁶ RAW 264.7 cells were intracellularly stained with 0.25 ug NOX2 Recombinant antibody (82872-1-RR, Clone:1B17) and CoraLite®488-Conjugated Goat Anti-Rabbit IgG(H+L) (SA00013-2)(red), or 0.25 ug Rabbit IgG Isotype Control Recombinant Antibody (98136-1-RR, Clone: 240953C9) (blue). Cells were fixed with 4% PFA and permeabilized with Flow Cytometry Perm Buffer (PF00011-C). This data was developed using the same antibody clone with 82872-1-PBS in a

Immunofluorescent analysis of (4% PFA) fixed paraffin-embedded mouse spleen tissue using NOX2 antibody (82872-1-RR, Clone: 1B17) at dilution of 1:2000 and CoraLite®488-Conjugated AffiniPure Goat Anti-Rabbit IgG(H+L) (SA00013-2). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 82872-1-PBS in a different storage buffer formulation.