

For Research Use Only

APOL6 Recombinant antibody

Catalog Number: 82877-1-RR



Basic Information

Catalog Number:

82877-1-RR

Size:

1000 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG24430

GenBank Accession Number:

BC038950

GeneID (NCBI):

80830

UNIPROT ID:

Q9BWW8

Full Name:

apolipoprotein L, 6

Calculated MW:

38 kDa

Purification Method:

Protein A purification

CloneNo.:

230162C6

Applications

Tested Applications:

ELISA

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

