

For Research Use Only

APOL6 Recombinant antibody

Catalog Number: 82877-1-RR



Basic Information

Catalog Number: 82877-1-RR	GenBank Accession Number: BC038950	Purification Method: Protein A purification
Size: 1000 µg/ml	GeneID (NCBI): 80830	CloneNo.: 230162C6
Source: Rabbit	UNIPROT ID: Q9BWW8	
Isotype: IgG	Full Name: apolipoprotein L, 6	
Immunogen Catalog Number: AG24430	Calculated MW: 38 kDa	

Applications

Tested Applications:
FC, ELISA

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1X10⁶ U2OS cells were intracellularly stained with 0.13 ug Anti-Human APOL6 (82877-1-RR, Clone:230162C6) and Alexa Fluor™ Plus 647-Conjugated AffiniPure Goat Anti-Rabbit IgG(H+L) at dilution 1:1000 (red), or 0.13 ug Rabbit IgG control Rabbit PolyAb (30000-0-AP, Clone:) (blue). Cells were fixed and permeabilized with True-Nuclear Transcription Factor Buffer Set.