

For Research Use Only

MMP8 Recombinant antibody, PBS Only (Detector)



Catalog Number: 82885-4-PBS

Basic Information

Catalog Number:

82885-4-PBS

Size:

1mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC074988

GeneID (NCBI):

4317

UNIPROT ID:

P22894

Full Name:

matrix metalloproteinase 8
(neutrophil collagenase)

Calculated MW:

467 aa, 53 kDa

Purification Method:

Protein A purification

CloneNo.:

230030G12

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve ofMP00013-1, MMP8 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82885-5-PBS. Detection antibody: 82885-4-PBS. Standard: Eg0378. Range: 0.156-10 ng/mL