

For Research Use Only

MMP8 Recombinant antibody, PBS Only (Capture/Detector)

Catalog Number: 82885-5-PBS



Basic Information

Catalog Number:

82885-5-PBS

Concentration:

1000 ug/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC074988

GeneID (NCBI):

4317

UNIPROT ID:

P22894

Full Name:

matrix metalloproteinase 8
(neutrophil collagenase)

Calculated MW:

467 aa, 53 kDa

Purification Method:

Protein A purification

CloneNo.:

230030H8

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA, Sample test

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Sandwich ELISA standard curve of MP00013-1, MMP8 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82885-5-PBS. Detection antibody: 82885-4-PBS. Standard: Eg0378. Range: 0.156-10 ng/mL



Sandwich ELISA standard curve of MP00013-2, MMP8 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82885-3-PBS. Detection antibody: 82885-5-PBS. Standard: Eg0378. Range: 0.156-10 ng/mL



Biolayer interferometry (BLI) kinetic assay of 82885-5-PBS against Human MMP8 was performed. The affinity constant is 3.31 nM.