

For Research Use Only

RBM17 Recombinant antibody

Catalog Number: 82889-3-RR



Basic Information

Catalog Number: 82889-3-RR	GenBank Accession Number: BC007871	Purification Method: Protein A purification
Size: 1000 µg/ml	GeneID (NCBI): 84991	CloneNo.: 230125A3
Source: Rabbit	UNIPROT ID: Q96125	Recommended Dilutions: WB 1:2000-1:10000
Isotype: IgG	Full Name: RNA binding motif protein 17	
Immunogen Catalog Number: AG0985	Calculated MW: 45 kDa	
	Observed MW: 45-50 kDa	

Applications

Tested Applications:
WB, ELISA

Species Specificity:
Human, Mouse

Positive Controls:

WB : BxPC-3 cells, A2780 cells, MCF-7 cells, HeLa cells, mouse thymus tissue

Background Information

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 82889-3-RR (RBM17 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours.