

For Research Use Only

TRIM62 Recombinant antibody, PBS Only (Capture)

Catalog Number: 82892-3-PBS



Basic Information

Catalog Number:

82892-3-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG7143

GenBank Accession Number:

BC001222

GeneID (NCBI):

55223

UNIPROT ID:

Q9BVG3

Full Name:

tripartite motif-containing 62

Calculated MW:

54 kDa

Purification Method:

Protein A purification

CloneNo.:

230195C11

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01312-2, TRIM62 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82892-3-PBS. Detection antibody: 82892-4-PBS. Standard: Ag7143. Range: 1.25-20 ng/mL.