

For Research Use Only

# FANCM Recombinant antibody

Catalog Number: 82897-1-RR



## Basic Information

Catalog Number:

82897-1-RR

Size:

1000  $\mu$ g/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG34447

GenBank Accession Number:

BC140776

GeneID (NCBI):

57697

UNIPROT ID:

Q8IYD8

Full Name:

Fanconi anemia, complementation group M

Purification Method:

Protein A purification

CloneNo.:

230233A12

Recommended Dilutions:

IF/ICC 1:750-1:3000

## Applications

Tested Applications:

IF/ICC, ELISA

Species Specificity:

human

Positive Controls:

IF/ICC : HeLa cells,

## Background Information

### Storage

Storage:

Store at  $-20^{\circ}\text{C}$ . Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for  $-20^{\circ}\text{C}$  storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Immunofluorescent analysis of (4% PFA) fixed HeLa cells using FANCM antibody (82897-1-RR, Clone: 230233A12) at dilution of 1:1500 and Multi-rAb CoraLite® Plus 488-Goat Anti-Rabbit Recombinant Secondary Antibody (H+L) (RGAR002), CL594-phalloidin (red).



Biolayer interferometry (BLI) kinetic assays of 82897-1-RR against Human FANCM were performed. The affinity constant is below 1 pM.