

For Research Use Only

# FANCM Recombinant antibody, PBS Only (Capture)

Catalog Number: 82897-3-PBS



## Basic Information

Catalog Number:

82897-3-PBS

Concentration:

1mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG34447

GenBank Accession Number:

BC140776

GeneID (NCBI):

57697

UNIPROT ID:

Q8IYD8

Full Name:

Fanconi anemia, complementation group M

Purification Method:

Protein A purification

CloneNo.:

230233A10

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP00057-3, FANCM Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82897-3-PBS. Detection antibody: 82897-4-PBS. Standard: Ag34447. Range: 0.313-20 ng/mL