

For Research Use Only

# FANCM Recombinant antibody, PBS Only (Detector)



Catalog Number: 82897-5-PBS

## Basic Information

Catalog Number:

82897-5-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG34447

GenBank Accession Number:

BC140776

GeneID (NCBI):

57697

UNIPROT ID:

Q8IYD8

Full Name:

Fanconi anemia, complementation

group M

Purification Method:

Protein A purification

CloneNo.:

230233D10

## Applications

Tested Applications:

Cytometric bead array, ELISA

Species Specificity:

Human

## Background Information

### Storage

Storage:

Store at -80°C.

Storage Buffer:

100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Standard curve of MP00057-1 measured by cytometric bead array. Capture antibody:82897-6-PBS. Detection antibody: 82897-5-PBS. Range: 0.313-20 ng/mL