

For Research Use Only

FANCM Recombinant antibody, PBS Only (Detector)



Catalog Number: 82897-5-PBS

Basic Information

Catalog Number:

82897-5-PBS

Size:

1mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG34447

GenBank Accession Number:

BC140776

GeneID (NCBI):

57697

UNIPROT ID:

Q8IYD8

Full Name:

Fanconi anemia, complementation group M

Purification Method:

Protein A purification

CloneNo.:

230233D10

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00057-1, FANCM Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82897-6-PBS. Detection antibody: 82897-5-PBS. Standard: Ag34447. Range: 0.313-20 ng/mL