

For Research Use Only

# HLA-DMA Recombinant antibody, PBS Only (Capture)



Catalog Number: 82899-5-PBS

## Basic Information

Catalog Number: 82899-5-PBS	GenBank Accession Number: BC026279	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 3108	CloneNo.: 230155F5
Source: Rabbit	UNIPROT ID: Q31604	
Isotype: IgG	Full Name: major histocompatibility complex, class II, DM alpha	
Immunogen Catalog Number: AG2452	Calculated MW: 261 aa, 29 kDa	

## Applications

Tested Applications:  
Indirect ELISA, Cytometric bead array

Species Specificity:  
Human

## Background Information

### Storage

Storage:  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:  
100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve ofMP00018-1, HLA-DMA Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82899-5-PBS. Detection antibody: 82899-1-PBS. Standard: Ag2452. Range: 0.625-20 ng/mL



Cytometric bead array standard curve ofMP00018-2, HLA-DMA Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82899-5-PBS. Detection antibody: 82899-4-PBS. Standard: Ag2452. Range: 0.625-20 ng/mL