

For Research Use Only

MET Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number: 82904-2-PBS



Basic Information

Catalog Number:

82904-2-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG0230

GenBank Accession Number:

BC130420

GeneID (NCBI):

4233

UNIPROT ID:

P08581

Full Name:

met proto-oncogene (hepatocyte
growth factor receptor)

Calculated MW:

1390 aa, 155 kDa

Purification Method:

Protein A purification

CloneNo.:

253142C9

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean MET concentration was determined to be 40.2 ng/mL in A431 cell extract based on a 1.5 mg/mL extract load, SH-SY5Y cell extract were served as a negative control.



Cytometric bead array standard curve of MP03178-1, MET Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82904-3-PBS. Detection antibody: 82904-2-PBS. Standard: Eg0230. Range: 0.156-20 ng/mL.