

For Research Use Only

BAZ2A Recombinant antibody, PBS Only (Capture)

Catalog Number: 82910-5-PBS



Basic Information

Catalog Number: 82910-5-PBS	GenBank Accession Number: NM_013449.4	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 11176	CloneNo.: 230205H2
Source: Rabbit	UNIPROT ID: Q9UIF9	
Isotype: IgG	Full Name: bromodomain adjacent to zinc finger domain, 2A	
Immunogen Catalog Number: AG33963	Calculated MW: 211 kDa	

Applications

Tested Applications:
Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01313-1, BAZ2A Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82910-5-PBS. Detection antibody: 82910-1-PBS. Standard: Ag33963. Range: 0.625-20 ng/mL