

For Research Use Only

TBCD Recombinant antibody, PBS Only

Catalog Number: 82911-2-PBS



Basic Information

Catalog Number:

82911-2-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG6310

GenBank Accession Number:

BC003094

GeneID (NCBI):

6904

UNIPROT ID:

Q9BTW9

Full Name:

tubulin folding cofactor D

Calculated MW:

133 kDa

Purification Method:

Protein A purification

CloneNo.:

230209D7

Applications

Tested Applications:

Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Biolayer interferometry (BLI) kinetic assays of 82911-2-RR against Human TBCD were performed. The affinity constant is 2.07 nM.