

For Research Use Only

# TBCD Recombinant antibody

Catalog Number: 82911-2-RR



## Basic Information

Catalog Number:

82911-2-RR

Size:

300 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG6310

GenBank Accession Number:

BC003094

GeneID (NCBI):

6904

UNIPROT ID:

Q9BTW9

Full Name:

tubulin folding cofactor D

Calculated MW:

133 kDa

Purification Method:

Protein A purification

CloneNo.:

230209D7

## Applications

Tested Applications:

ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Biolayer interferometry (BLI) kinetic assays of 82911-2-RR against Human TBCD were performed. The affinity constant is 2.07 nM.