

For Research Use Only

# CTNNBIP1 Recombinant antibody

Catalog Number: 82914-1-RR



## Basic Information

Catalog Number:	82914-1-RR	GenBank Accession Number:	BC014300	Purification Method:	Protein A purification
Size:	500 µg/ml	GeneID (NCBI):	56998	CloneNo.:	230164H5
Source:	Rabbit	UNIPROT ID:	Q9NSA3		
Isotype:	IgG	Full Name:	catenin, beta interacting protein 1		
Immunogen Catalog Number:	AG13160	Calculated MW:	21 kDa		

## Applications

Tested Applications:  
ELISA

Species Specificity:  
Human

## Background Information

### Storage

Storage:  
Store at -20°C. Stable for one year after shipment.

Storage Buffer:  
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Biolayer interferometry (BLI) kinetic assays of 82914-1-RR against Human CTNNBIP1 were performed. The affinity constant is 15.6 nM.