

For Research Use Only

Alpha-1-Antitrypsin Recombinant antibody, PBS Only (Detector)

Catalog Number: 82918-2-PBS



Basic Information

Catalog Number:

82918-2-PBS

Size:

1mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG9369

GenBank Accession Number:

BC015642

GeneID (NCBI):

5265

UNIPROT ID:

P01009

Full Name:

serpin peptidase inhibitor, clade A
(alpha-1 antiproteinase, antitrypsin),
member 1

Calculated MW:

418 aa, 47 kDa

Purification Method:

Protein A purification

CloneNo.:

230109C11

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA, Sample test

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP00022-1, Alpha-1-Antitrypsin Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82918-3-PBS. Detection antibody: 82918-2-PBS. Standard: Ag9369. Range: 0.39-25 ng/mL.