

For Research Use Only

NFX1 Recombinant antibody, PBS Only

Catalog Number: 82924-1-PBS



Basic Information

Catalog Number:

82924-1-PBS

Size:

1mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG9448

GenBank Accession Number:

BC012151

GeneID (NCBI):

4799

UNIPROT ID:

Q12986

Full Name:

nuclear transcription factor, X-box binding 1

Calculated MW:

1120 aa, 124 kDa

Purification Method:

Protein A purification

CloneNo.:

230211D9

Applications

Tested Applications:

IF/ICC, FC (Intra), Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:



T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1x10⁶ PC-3 cells were intracellularly stained with 0.25 ug NFX1 Recombinant antibody (82924-1-RR, Clone:230211D9) and APC-Conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)(red), or 0.25 ug Isotype Control (blue). Cells were fixed and permeabilized with True-Nuclear Transcription Factor Buffer Set. This data was developed using the same antibody clone with 82924-1-PBS in a different storage buffer formulation.

Immunofluorescent analysis of (4% PFA) fixed HeLa cells using NFX1 antibody (82924-1-RR, Clone: 230211D9) at dilution of 1:400 and CoraLite® 488-Conjugated AffiniPure Goat Anti-Rabbit IgG(H+L) (SA00013-2). This data was developed using the same antibody clone with 82924-1-PBS in a different storage buffer formulation.