

For Research Use Only

# RORA Recombinant antibody, PBS Only (Capture)

Catalog Number: 82930-6-PBS



## Basic Information

Catalog Number:

82930-6-PBS

Concentration:

1mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG0944

GenBank Accession Number:

BC008831

GeneID (NCBI):

6095

UNIPROT ID:

P35398

Full Name:

RAR-related orphan receptor A

Calculated MW:

59 kDa

Purification Method:

Protein A purification

CloneNo.:

230217E4

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP00055-2, RORA Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82930-6-PBS. Detection antibody: 82930-2-PBS. Standard: Ag0944. Range: 0.625-20 ng/mL



Cytometric bead array standard curve of MP00055-3, RORA Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82930-6-PBS. Detection antibody: 82930-4-PBS. Standard: Ag0944. Range: 0.625-20 ng/mL