

For Research Use Only

TAT Recombinant antibody

Catalog Number: 82939-1-RR



Basic Information

Catalog Number:	82939-1-RR	GenBank Accession Number:	BC130534	Purification Method:	Protein A purification
Size:	1000 µg/ml	GeneID (NCBI):	6898	CloneNo.:	230314H2
Source:	Rabbit	UNIPROT ID:	P17735		
Isotype:	IgG	Full Name:	tyrosine aminotransferase		
Immunogen Catalog Number:	AG34043	Observed MW:	50-60 kDa		

Applications

Tested Applications:
ELISA

Species Specificity:
Human, mouse

Background Information

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

