

For Research Use Only

# TAT Recombinant antibody, PBS Only (Capture)

Catalog Number: 82939-3-PBS



## Basic Information

Catalog Number:

82939-3-PBS

Size:

1mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG34043

GenBank Accession Number:

BC130534

GeneID (NCBI):

6898

UNIPROT ID:

P17735

Full Name:

tyrosine aminotransferase

Purification Method:

Protein A purification

CloneNo.:

230314F7

## Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Sandwich ELISA standard curve of MP00131-1, TAT Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82939-3-PBS. Detection antibody: 82939-2-PBS. Standard: Ag34043. Range: 0.313-20 ng/mL



Biolayer interferometry (BLI) kinetic assay of 82939-3-PBS against Human TAT was performed. The affinity constant is 0.981 nM.