

For Research Use Only

CCR5 Recombinant antibody, PBS Only



Catalog Number: 82942-1-PBS

Basic Information

Catalog Number: 82942-1-PBS	GenBank Accession Number: BC038398	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 1234	CloneNo.: 230266A6
Source: Rabbit	UNIPROT ID: P51681	
Isotype: IgG	Full Name: chemokine (C-C motif) receptor 5	
Immunogen Catalog Number: AG11493	Calculated MW: 352 aa, 41 kDa	
	Observed MW: 45 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
Human, mouse

Background Information

CCR5 (also known as CD195) is a member of the beta chemokine receptor family. This protein is expressed by T cells and macrophages and is known to be an important co-receptor for macrophage-tropic viruses, including HIV, to enter host cells. The natural chemokine ligands that bind to CCR5 include MIP-1-alpha, MIP-1-beta, MCP-2, and RANTES.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 82942-1-RR (CCR5 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 82942-1-PBS in a different storage buffer formulation.



Biolayer interferometry (BLI) kinetic assays of 82942-1-RR against Human CCR5 were performed. The affinity constant is 0.921 nM.