

For Research Use Only

PDGFR beta/CD140b Recombinant antibody, PBS Only (Detector)

Catalog Number: 82943-3-PBS



Basic Information

Catalog Number:

82943-3-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC032224

GeneID (NCBI):

5159

Full Name:

platelet-derived growth factor receptor, beta polypeptide

Calculated MW:

1106 aa, 124 kDa

Purification Method:

Protein A purification

CloneNo.:

230214E7

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA, Sample test

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP00006-3, PDGFR beta/CD140b Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82943-6-PBS. Detection antibody: 82943-3-PBS. Standard: Eg0147. Range: 0.156-10 ng/mL



Sandwich ELISA standard curve of MP00006-4, PDGFR beta Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82943-7-PBS. Detection antibody: 82943-3-PBS. Standard: Eg0147. Range: 0.156-10 ng/mL