

For Research Use Only

IFT140 Recombinant antibody, PBS Only (Detector)

Catalog Number: 82947-6-PBS



Basic Information

Catalog Number:

82947-6-PBS

Concentration:

1mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG11119

GenBank Accession Number:

BC035577

GeneID (NCBI):

9742

UNIPROT ID:

Q96RY7

Full Name:

intraflagellar transport 140 homolog
(Chlamydomonas)

Calculated MW:

1462 aa, 165 kDa

Purification Method:

Protein A purification

CloneNo.:

230231H2

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00029-2, IFT140 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82947-4-PBS. Detection antibody: 82947-6-PBS. Standard: Ag11119. Range: 2.5-20 ng/mL.



Cytometric bead array standard curve of MP00029-3, IFT140 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82947-3-PBS. Detection antibody: 82947-6-PBS. Standard: Ag11119. Range: 2.5-20 ng/mL.