

For Research Use Only

# TAB2 Recombinant antibody, PBS Only (Capture)

Catalog Number: 82948-3-PBS



## Basic Information

Catalog Number:

82948-3-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG5742

GenBank Accession Number:

BC035910

GeneID (NCBI):

23118

UNIPROT ID:

Q9NYJ8

Full Name:

mitogen-activated protein kinase  
kinase kinase 7 interacting protein 2

Calculated MW:

76 kDa

Purification Method:

Protein A purification

CloneNo.:

230118C6

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP00082-1, TAB2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82948-3-PBS. Detection antibody: 82948-2-PBS. Standard: Ag5742. Range: 0.625-20 ng/mL



Cytometric bead array standard curve of MP00082-3, TAB2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82948-3-PBS. Detection antibody: 82948-1-PBS. Standard: Ag5742. Range: 0.625-20 ng/mL