

For Research Use Only

GPX4 Recombinant antibody, PBS Only (Detector)

Catalog Number: 82949-2-PBS

1 Publications



Basic Information

Catalog Number:

82949-2-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG5809

GenBank Accession Number:

BC022071

GeneID (NCBI):

2879

UNIPROT ID:

P36969

Full Name:

glutathione peroxidase 4
(phospholipid hydroperoxidase)

Calculated MW:

22 kDa

Purification Method:

Protein A purification

CloneNo.:

230168B5

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Cited Applications:

WB, IF

Species Specificity:

human

Background Information

Notable Publications

Author	Pubmed ID	Journal	Application
Lin Luo	39667244	Tissue Cell	WB,IF

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00049-2, GPX4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82949-3-PBS. Detection antibody: 82949-2-PBS. Standard: Ag5809. Range: 2.5-20 ng/mL.