

For Research Use Only

EFNA2 Recombinant antibody, PBS Only (Capture)



Catalog Number: 82951-5-PBS

Basic Information

Catalog Number:

82951-5-PBS

Size:

1mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_001405

GeneID (NCBI):

1943

UNIPROT ID:

O43921

Full Name:

ephrin-A2

Calculated MW:

213 aa, 24 kDa

Purification Method:

Protein A purification

CloneNo.:

230251H4

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve ofMP00081-1, EFNA2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82951-5-PBS. Detection antibody: 82951-4-PBS. Standard: SY00429. Range: 0.313-10 ng/mL



Cytometric bead array standard curve ofMP00081-2, EFNA2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82951-5-PBS. Detection antibody: 82951-3-PBS. Standard: SY00429. Range: 0.313-10 ng/mL