

For Research Use Only

uPAR/CD87 Recombinant antibody, PBS Only (Capture)

Catalog Number: 82953-4-PBS



Basic Information

Catalog Number:

82953-4-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG28056

GenBank Accession Number:

BC002788

GeneID (NCBI):

5329

UNIPROT ID:

Q03405

Full Name:

PLAUR

Calculated MW:

37 kDa

Purification Method:

Protein A purification

CloneNo.:

230148F8

Applications

Tested Applications:

Cytometric bead array, Sandwich ELISA, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP00129-2, uPAR/CD87 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82953-4-PBS. Detection antibody: 82953-2-PBS. Standard: Ag28056. Range: 125-8000 pg/mL



Cytometric bead array standard curve of MP00129-2, uPAR/CD87 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82953-4-PBS. Detection antibody: 82953-2-PBS. Standard: Ag28056. Range: 1.25-20 ng/mL