

For Research Use Only

BBS4 Recombinant antibody, PBS Only (Detector)

Catalog Number: 82961-3-PBS



Basic Information

Catalog Number:

82961-3-PBS

Concentration:

1mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG3430

GenBank Accession Number:

BC027624

GeneID (NCBI):

585

UNIPROT ID:

Q96RK4

Full Name:

Bardet-Biedl syndrome 4

Calculated MW:

519 aa, 58 kDa

Purification Method:

Protein A purification

CloneNo.:

230303E4

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Cytometric bead array standard curve of MP00130-1, BBS4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82961-1-PBS. Detection antibody: 82961-3-PBS. Standard: Ag3430. Range: 0.313-20 ng/mL.



Cytometric bead array standard curve of MP00130-2, BBS4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82961-2-PBS. Detection antibody: 82961-3-PBS. Standard: Ag3430. Range: 0.313-20 ng/mL.



Biolayer interferometry (BLI) kinetic assays of 82961-3-RR against Human BBS4 were performed. The affinity constant is 0.575 nM.