

For Research Use Only

# ACAP1 Recombinant antibody, PBS Only (Capture)

Catalog Number: 82967-2-PBS



## Basic Information

Catalog Number:

82967-2-PBS

Size:

1mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG34309

GenBank Accession Number:

BC018543

GeneID (NCBI):

9744

UNIPROT ID:

Q15027

Full Name:

ArfGAP with coiled-coil, ankyrin repeat and PH domains 1

Calculated MW:

740 aa, 82 kDa

Purification Method:

Protein A purification

CloneNo.:

230370G6

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP01212-1, ACAP1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82967-2-PBS. Detection antibody: 82967-3-PBS. Standard: Ag34309. Range: 1.25-20 ng/mL.